



SM-500D



15 inch color high resolution TFT display

Portability, Suitable for ICU, operating room monitor

Unique blood pressure measurement technique

Anti high-frequency electroartifacts and defibrillation, digital filter

Built-in rechargeable battery, DC/AC



SM-500D

15inchTFT



Picture:

- Big font display interface
- Trend graph display interface
- OxyCRG display interface
- ECG, NIBP, SPO2, TEMP interface
- Standard display interface
- Knob operate
- ECG full lead display interface

Feature:

- Display with 11 waveform
- 720 hours trend data
- 1000 NIBP data storage
- 200 alarm event data storage
- Sound, light, human voice alarm
- 12 hours waveform review
- Wire and wireless network
- Nurse call
- analog output
- Support software upgrade online
- Arrhythmia analysis, ST segment analysis
- Support touch screen

Standard:

ECG, RESP, NIBP, SpO2, TEMP

Optional:

Recorder, IBP, EtCO2, CSM, Multi-Gas, touch screen



Sonostar Technologies Co., Limited

Address: #31, Shuiwan road, Jida, Zhuhai, Guangdong, China
Tel: +86 13168355297 MSN: sonostar@163.com Skype: sonostar-ut
Web: www.sonostar.net Mail: sonostar@sonostar.net

