



*Dedicated to the popularization of ultrasonic technology*

## **SS-1** (for human)/**SS-1V** (for veterinary) **Palm Ultrasound Diagnosis B Scanner**

SS-1/SS-1V is enjoying a great success in global market by its innovation in technology and user-friendly design. It is distinguishable by the lightweight, easy hand carrying, accurate imaging and simple interface. It is really the most ideal for doctors, vet, even the farmer who can easily carry out quick and easy diagnosis for Uterine, Ovarian examination, Pregnancy detection, sexing and so on.



### **SYSTEM FEATURES:**

- Excellent image quality
- 5.0 inch TFT LCD
- Easy to operation
- Low power consumption
- Waterproof Probe

### **CLINICAL APPLICATIONS:**

Abdominal (human/veterinary)  
Obstetrical (human/veterinary)  
Wide veterinary application: swine, sheep, equine, cat, bovine, dog, etc.



#### PHYSICAL CHARACTERISTICS:

Dimension: 235\*130\*43 mm

Weight: 750g (excluding accessories)

Monitor	5.0 inch TFT LCD
Weight	750g(excluding accessories)
Probes	3.5MHz Mechanical Sector Selectable
Battery	Internal Rechargeable keeps working for 2h continually,4h for stand-by
Image Mode	B,B+B,B+M,M
Image Magnification	Image can be zoom in 4 steps:×1.0, ×1.2, ×1.6, ×2.0
Gray Scale	256 levels
Text Annotation	Patient ID, Sex, Age
Software	For human or veterinary, selectable
Measurement	Distance, Length, Area/Circumference, Time/Heart Rate/EF rate for heart
Stander Accessories	AC Adapter, Internal Battery, Gel, Mouse

#### STANDARD ACCESSORIES:

Main unit 1set

3.5MHz Mechanical Sector Selectable

Stander Accessories Adapter

Internal Battery

Pseudo color

Battery unit 1 set

Metallic case

#### OPTIONS:

Extra rechargeable internal battery with the cable mouse



## Sonostar Technologies Co., Limited

Address: #31, Shuiwan road, Jida, Zhuhai, Guangdong, China

Tel: 18925122158 MSN: sonostar@163.com Skype: sonostar-ut

Web: www.sonostar.net Mail: sonostar@sonostar.net

