

V3EC2014/A0



More information, please browse www.sonostar.net

Sonostar Technologies Co., Limited

Headquarters: C Building, #27 Yayingshi Road, Science Town, Guangzhou, China
 Parts Factory: First floor, Investment Building, Nanping, Zhuhai, Guangdong
 Tel: +86-20-32382095 Fax: +86-20-62614030 Skype: sonostar-ut
www.sonostar.net Mail: sonostar@sonostar.net MSN: sonostar@163.com

Note: Sonostar trademarks is protected by Madrid Treaty.



Stable & Good Image!

V3

Palm Ultrasound Scanner

- 5.6 inch LCD screen with high resolution image
- Stable case
- Simple and easy-touch keyboard
- Built-in Lithium battery, long lasting
- Large capacity SD card for memory
- Waterproof, easy to clean





Rectal probe for cattle & horse



Specification:

- Screen: 5.6 inches
- Display mode: B
- Scan: Sector Scan
- Probe: 3.5MHz/7MHz optional
- Scan Depth: 95-185mm
- Gain: 30dB-105dB
- Image Store: 8GB (SD Card)
- Gray Scale: 256
- Battery last: 8 hours
- Size: 210mmX120mmX40mm
- Weight: 1.2kg with battery
- Probe Connector: 1

Standard Configuration:

- Main Host 1 Unit
- 3.5MHz/5.0MHz Mechanical Sector Selectable 1 Pcs
- Battery charger 1 Set
- SD card 1 Pcs

Optional:

- Rectal Pobe
- Battery
- Car charger



Waterproof